



GRAHA TEKNO MEDIKA

Hospital Furniture & Equipment Manufacturing



EMERGENCY TROLLEY 3 DRAWER 8211801



SPECIFICATIONS :

- | | |
|---------------------|---|
| Construction | : Frame made of steel.
Consists of 3 drawers and 1 cabinet.
Handle made of stainless steel. |
| Dimension | : 700 (L) x 450 (W) x 900 (H) mm. |
| Castor | : Four 4" castors, two with brakes. |
| Finishing | : Powder coating. |
| Accessories | : Folding table made of steel.
IV stand made of stainless steel.
Oxygen holder made of steel. |



Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com

www.grahateknomedika.com

MADE IN INDONESIA

