



LINEN HAMPER CARRIAGE 8211702



SPECIFICATIONS :

- Construction** : Frame made of steel pipe and plate.
Consists of three shelves.
Handle made of steel.
- Dimension** : 100 (L) x 65 (W) x 110 (H) mm.
- Castors** : Four 4" castors, two with brakes.
- Finishing** : Powder coating.
- Accessories** : Hamper made of canvas.



MADE IN INDONESIA

Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com

www.grahateknomedika.com