



BED SCREEN **2 SCREEN** 8210602

SPECIFICATIONS :

Construction :

Frame made of round pipe and box pipe.
Panel consists of two pieces.
Screen made from cotton.

Dimension :

1822 (L) x 300 (W) x 1700 (H) mm.

Castors :

Five 2" castors.

Finishing :

Powder Coating.



Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com

www.grahateknomedika.com

MADE IN INDONESIA

