



Gris
HOSPITAL EQUIPMENT

DONOR TABLE 8211010

SPECIFICATIONS :

- | | |
|-------------------------|--|
| Construction | : Main frame made of steel pipe.
Back rest frame made of steel pipe. |
| Dimension | : 2000 (L) x 600 (W) x 800 (H) mm. |
| Adjustment Range | : Adjustable by gas spring up to 70°. |
| Castor | : Four 4" castors, two with brakes. |
| Finishing | : Powder coating. |
| Accessories | : Mattress made of foam covered by vinyl.
Hand support 300 x 200 mm with steel elbow. |



Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com
www.grahateknomedika.com

MADE IN INDONESIA

