



Gris
HOSPITAL EQUIPMENT

DRESSING TROLLEY

1 DRAWER
8220502

SPECIFICATIONS :

- Construction** : Frame made of stainless steel.
Handle made of stainless steel.
Drawer made of stainless steel.
- Dimension** : 600 (L) x 400 (W) x 800 (H) mm.
- Castor** : Four 3" castors, two with brakes.
- Finishing** : Polished.
- Accessories** : Bowl made of stainless steel.
Bowl diameter 30 cm.



Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com

www.grahateknomedika.com

MADE IN INDONESIA

