



EMERGENCY TROLLEY

STAINLESS STEEL
3 DRAWER
8221801



SPECIFICATIONS :

- Construction** : Frame made of stainless steel.
Consists of 3 drawers and 1 cabinet.
Handle made of stainless steel.
- Dimension** : 700 (L) x 450 (W) x 900 (H) mm.
- Castor** : Four 4" castors, two with brakes.
- Finishing** : Polished.
- Accessories** : Folding table made of stainless steel.
IV stand made of stainless steel.
Oxygen holder made of stainless steel.



MADE IN INDONESIA

Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com

www.grahateknomedika.com