



EMERGENCY STRETCHER BACK RAISE

8211004



SPECIFICATIONS :

- | | |
|-------------------------|--|
| Construction | : Upper frame made of steel pipe.
Back rest frame made of steel pipe.
Side guard made of steel pipe. |
| Dimension | : 2000 (L) x 600 (W) x 800 (H) mm. |
| Adjustment Range | : Back rest adjustable by gas spring up to 70°. |
| Castors | : Four 4" castors, two with brakes. |
| Finishing | : Powder coating. |
| Accessories | : Oxygen holder made of steel pipe.
IV stand made of stainless steel.
Mattress made of foam covered by vinyl.
Multi-purpose tray. |



Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com

www.grahateknomedika.com

MADE IN INDONESIA

