



EXAMINATION TABLE STAINLESS STEEL BACK RAISE 8221006



SPECIFICATIONS :

- | | |
|-------------------------|--|
| Construction | : Frame made of stainless steel pipe.
Back rest frame made of stainless steel pipe. |
| Dimension | : 2000 (L) x 600 (W) x 800 (H) mm. |
| Adjustment Range | : Back rest adjustable by gas spring up to 70°. |
| Finishing | : Polished. |
| Accessories | : Mattress made of foam covered by vinyl.
Foot step made of stainless steel.
Foot step dimension 440 (L) x 440 (W) x 410 (H) mm. |



Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com

www.grahateknomedika.com

MADE IN INDONESIA

