




PRODUCT KEY BENEFITS

- ABS plastic headboard and footboard with modern and sleek design.
- Back rest retraction to ensure patient's comfort and ergonomic position.
- Standard side rail to minimize risk of falling, easy to clean and to disinfect.
- Columns and actuators to ensure good and reliable performance.
- Two central locks at head side and foot side, to ensure safety maneuver while moving the bed.
- Hi - low function.
- Adjustable trendelenburg and reverse trendelenburg position.
- Adjustable tilting position.


SPECIFICATIONS OF ICU BED 5 AXIS AB0115

Overall Dimensions	: 2160 mm x 1085 mm x 1005 mm.
Construction	: Made of 30 mm x 60 mm and 30 mm x 30 mm steel pipe.
Positions	: Operated with electric motor system. Include remote control (220 VAC - 50 Hz) Ergonomic back rest design with retraction up to 110 mm and elevation up to 75° Adjustable knee rest with elevation up to 35° Adjustable trendelenburg elevation up to 18°, -18° Adjustable tilting elevation up to 13°, -13° Adjustable high low position with range of 605 mm - 1005 mm without mattress CPR position.
Mattress Platform	: Made of steel sheet with hollow steel pipe for safety frame.
Head & Foot Panel	: Made of ABS plastic with bumper, detachable.
Safety Side Rail	: Made of ABS plastic, collapsible. Automatic locking and unlocking system with one touch folding.
Castor	: Four 5" castors, with two central locking system.
IV Stand	: Made of stainless steel pipe (with two hanger hook), 6 optional locations.
Painting	: Powder Coating.
Accessories	: Urinal bag hanger (stainless steel), Waffle mattress (2000 x 905 x 120 mm).
Maximum Safe Working Load	: 250 kgs.
Optional Features	: Weighing Scale, Standard Head & Foot Panel Model 2, Standard Safety Side Rail Model 2.

Standard Features :



Optional Features :



Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com

www.grahateknomedika.com