



GRAHA TEKNO MEDIKA

Hospital Furniture & Equipment Manufacturing



BEDSIDE CABINET

8210307




SPECIFICATIONS :

Construction : Frame made of steel plate.
Upper table made of HPL board.
One door, made of HPL board completed with lock and handle.
One drawer, made of HPL board completed with lock and handle.
Removeable shoe shelf part.
Stainless Steel towel hanger.

Finishing : Powder coating.

Castor : 2".

Dimension : 450 (L) x 400 (W) x 830 (H) mm.



Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com

www.grahateknomedika.com

MADE IN INDONESIA

