



GRAHA TEKNO MEDIKA

Hospital Furniture & Equipment Manufacturing



WAFFLE MATTRESS

999002

SPECIFICATIONS :

- Mattress Type** : Waffle mattress.
Density : 23.
Cover : Local anti microbial PVC.
Dimension : 2000 (L) x 840 (W) x 120 (H) mm.



MADE IN INDONESIA

Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com

www.grahateknomedika.com