

SPECIFICATIONS OF INFANT WARMER - 921212

Caregiver Controls	: 4" monochrome touch screen panel.
Warming Mode	: 27°C - 38°C, 1°C resolution.
Air Mode	: 1°C.
Display Resolution	: Less than 30 minutes.
Warming Up Time	: Temperature out of limit.
Safety Alarms	Sensor defect. Sensor detached.
Additional Alarms	: Appar time 1,3,5,10 minutes.
Safety (Automatic Off)	: Air more than 40°C. Temperature setup. Sensor failure. 3°C deviation from setting.
Examination Lamp	: 3 watt LED.
Construction Material (Frame)	: Steel sheet and steel tubing or powder coating.
Mattress	
Mattress size	: 73 x 53 x 5 cm.
Mattress height from ground	: Approx 95 cm.
Mattress slope	: Adjustable manually from -12° to 12°.
IV Pole	: Adjustable with fixation clamp.
Heater	: 500 watt quartz heater.
Temperature Controller	: PID System.
Temperature Sensor	: Air Sensor.
Overall Dimensions	: 110 (L) x 60 (W) x 190 (H) cm.
Castor	
Diameter	: 5".
Swivel Locks	: 360°.
Brakes	: Castor with stopper.
Brakes type	: Two castors with lock.
Max Patient Weight	: 10 Kg.
Electrical	
Power Required, VAC	: 220 V , 50 Hz.
Purchase Information	
Warranty	: 1 Year.
Spare Part	: 5 Years.

Features :



Bassinet adjustable to Trendelenburg and reverse.



IV Stand

PRODUCT KEY BENEFITS

- Programmable logic controller with full touch screen.
- Touch screen display can be customized according to customer preference.
- High level of safety performance.
- Input bias mode to adjust the warmer against unsuitable environments.



4" monochrome touch screen controller

Marketing Office :
 Jl. Sultan Iskandar Muda Kav. 24
 Jakarta Selatan 12240 - Indonesia
 P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :
 Jl. Abdul Wahab No. 59, Cinangka,
 Sawangan, Depok - Indonesia
 P. +6221 7494512 | F. +6221 7494512
info@grahateknomedika.com
www.grahateknomedika.com