




INFUSION STAND 8211201

SPECIFICATIONS :

- | | |
|-------------------------|--|
| Construction | : Upper stand made of stainless steel pipe.
Lower stand made of steel pipe.
Two hooks made of stainless steel. |
| Dimension | : 600 x 1200 - 2200 (H) mm. |
| Adjustment Range | : Hi-Low adjustment by screw mechanism. |
| Castor | : Four 2" castors. |
| Finishing | : Powder coating and polished. |



MADE IN INDONESIA

Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com
www.grahateknomedika.com