



INFUSION STAND

STAINLESS STEEL

8221201

SPECIFICATIONS :

- Construction** : Upper stand made of stainless steel pipe.
Lower stand made of stainless steel pipe.
Two hooks made of stainless steel.
- Dimension** : 600 x 1200 - 2200 (H) mm.
- Adjustment Range** : Hi-Low adjustment by screw mechanism.
- Castor** : Four 2" castors.
- Finishing** : Polished.



Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com
www.grahateknomedika.com

MADE IN INDONESIA

