



**INSTRUMENT
CABINET**
STAINLESS STEEL
2 DOOR FULL GLASS
8221302



SPECIFICATIONS :

- Construction** : Main frame made of stainless steel pipe.
Upper frame made of stainless steel plate.
Four shelves made of glass (5 mm).
Door made of stainless steel square pipe (completed with lock and two handles).
- Dimension** : 1200 (L) x 450 (W) x 1700 (H) mm.
- Castors** : Four 3" castors, two with brakes.
- Finishing** : Polished.

Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com
www.grahateknomedika.com

MADE IN INDONESIA

