



MAYO TABLE

STAINLESS STEEL
8220902



SPECIFICATIONS :

Construction	: Table made of stainless steel. Upper stand made of stainless steel. Lower stand made of stainless steel.
Dimension	: 700 (L) x 400 (W) x 740 - 1150 (H) mm.
Adjustment Range	: Hi-Low adjustment by screw mechanism.
Castor	: Four 2" castors.
Finishing	: Polished



Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com

www.grahateknomedika.com

MADE IN INDONESIA

