



GRAHA TEKNO MEDIKA

Hospital Furniture & Equipment Manufacturing



OVERBED TABLE 8210901

SPECIFICATIONS :

Contruction	: Upper stand made of stainless steel.
	: Lower stand made of stainless steel.
	: Foot unit made of steel pipes.
Table	: Plywood covered by melamine.
Adjustment Range	: Hi low adjustment by screw mecanism.
Finishing	: Powder coating and polished.
Castor	: 2".
Dimension	: 700 (L) x 400 (w) x 740 - 1150 (H) mm.



Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

Info@grahateknomedika.com
www.grahateknomedika.com



MADE IN INDONESIA