



PHOTOTHERAPY

8212101



SPECIFICATIONS :

Contruction	: Main frame made of square steel plate. : Cover made of steel plate. : Foot unit made of square steel plate.
Lamp	: 4 pcs x 20 watt blue light lamp.
Time Control	: Using hour parameter.
Finishing	: Powder coating.
Castor	: 2".
Dimension	: 690 (L) x 690 (W) x 1320 - 2000 (H) mm.



Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com
www.grahateknomedika.com

MADE IN INDONESIA

