



GRAHA TEKNO MEDIKA

Hospital Furniture & Equipment Manufacturing



BABY BOX (ACRYLIC)

8220701



SPECIFICATIONS :

- | | |
|------------------------|---|
| Contruction | : Basket made of stainless steel rod and acrylic. |
| | : Trolley made of stainless steel pipes. |
| Finishing | : Polished. |
| Castor | : 2". |
| Dimension | : Basket - 630 x 390 x 370 mm. |
| (L) x (W) x (H) | : Trolley - 700 x 440 x 800 mm. |
| Accessories | : Mattress made of foam covered by vinyl. |



Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com

www.grahateknomedika.com

MADE IN INDONESIA

