



BEDSIDE CABINET 8210301



SPECIFICATIONS :

- | | |
|--------------------|---|
| Contruction | : Frame made of steel plate. |
| | : Upper table made of steel plate. |
| | : One door, made of steel plate completed with lock and handle. |
| | : One drawer, made of steel plate completed with lock and handle. |
| Finishing | : Powder coating. |
| Castor | : 2". |
| Dimension | : 600 (L) x 400 (W) x 800 (H) mm. |



Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com
www.grahateknomedika.com

MADE IN INDONESIA

