



PRODUCT KEY BENEFITS

- HPL head board and foot board are easier to clean and to disinfect.
- HPL color can be customized according to the room interior.
- Back rest retraction to ensure patient's comfort and ergonomic position.
- Stainless side rail to minimize risk of falling, no entrapment zone, easy to clean and to disinfect.
- Foldable crank with torque release system.
- Upgradeable to manual 2 crank, 3 crank, and electric hospital bed.

SPECIFICATIONS OF HOSPITAL BED 1 CRANK - 0121

Overall Dimensions	: 2160 mm x 980 mm x 635 mm.
Construction	: Made of 30 mm x 60 mm and 30 mm x 30 mm steel pipe.
Crank Mechanism	: Ergonomic back rest design with retraction up to 110 mm and elevation up to 75°.
Mattress Platform	: Made of steel sheet with hollow steel pipe safety frame.
Head & Foot Panel	: Made of High Pressure Laminate (HPL), detachable.
Safety Side Rail	: Made of stainless steel, collapsible. Automatic locking and unlocking system with one touch folding.
Castor	: Four 5" castors, 360° swivel, with 2 brakes.
IV Stand	: Made of stainless steel pipe (with two hanger hooks), 6 optional locations.
Painting	: Powder coating.
Accessories	: Name tag (acrylic), Urinal bag hanger (stainless steel), Mattress (superfoam 2000 x 840 x 120 mm).
Maximum Safe Working Load	: 250 Kgs.

Optional :

1. Upgradeable to Hospital Bed 1 Axis Electric.
2. Upgradeable to Hospital Bed 2 Crank and 2 Axis Electric.
3. Upgradeable to Hospital Bed 3 Crank and 3 Axis Electric (by replacing the lower deck).

Standard Features :



IV Stand




Urinal Bag
Hanger



European
Standard Side Rail



Wide Name
Tag Holder



Standard
Mattress

Color Option :



Wooden



Blue



Yellow



Green

Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com

www.grahateknomedika.com