



## PRODUCT KEY BENEFITS

- ABS plastic head board and foot board with modern and sleek design.
- Back rest retraction to ensure patient's comfort and ergonomic position.
- Standard side rail to minimize risk of falling, easy to clean and to disinfect.
- Foldable crank with torque release system.
- Upgradeable to manual 2 crank, 3 crank and electric hospital bed.

## SPECIFICATIONS OF HOSPITAL BED 1 CRANK ABS - SB 0124

|                                  |                                                                                                           |
|----------------------------------|-----------------------------------------------------------------------------------------------------------|
| <b>Overall Dimensions</b>        | : 2160 mm x 980 mm x 635 mm.                                                                              |
| <b>Construction</b>              | : Made of 30 mm x 60 mm and 30 mm x 30 mm steel pipe.                                                     |
| <b>Crank Mechanism</b>           | : Ergonomic back rest design with retraction up to 110 mm and elevation up to 75°.                        |
| <b>Mattress Platform</b>         | : Made of steel sheet with hollow steel pipe safety frame.                                                |
| <b>Head &amp; Foot Panel</b>     | : Made of ABS plastic with bumper, detachable.                                                            |
| <b>Safety Side Rail</b>          | : Made of aluminium alloy, collapsible.<br>Automatic locking and unlocking system with one touch folding. |
| <b>Castor</b>                    | : Four 5" castors, 360° swivel, with 2 brakes.                                                            |
| <b>IV Stand</b>                  | : Made of stainless steel pipe (with two hanger hooks), 6 optional locations.                             |
| <b>Painting</b>                  | : Powder coating.                                                                                         |
| <b>Accessories</b>               | : Urinal bag hanger (stainless steel), Mattress (superfoam 2000 x 840 x 120 mm).                          |
| <b>Maximum Safe Working Load</b> | : 250 Kgs.                                                                                                |

### Optional :

1. Upgradeable to Hospital Bed 1 Axis Electric.
2. Upgradeable to Hospital Bed 2 Crank and 2 Axis Electric.
3. Upgradeable to Hospital Bed 3 Crank and 3 Axis Electric (by replacing the lower deck).

### Standard Features :



IV Stand



Urinal Bag Hanger



Standard Side Rail



Standard Mattress

### Color Option :



Blue



Green



Yellow



Wooden

### Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24  
 Jakarta Selatan 12240 - Indonesia  
 P. +6221 7264150 | F. +6221 7264152

### Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,  
 Sawangan, Depok - Indonesia  
 P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com

www.grahateknomedika.com