



INSTRUMENT CABINET

2 DOOR FULL GLASS
8211302



SPECIFICATIONS :

- Construction** : Main frame made of steel pipe.
Upper frame made of steel plate.
Four shelves made of glass (5 mm).
Door made of square steel pipe,
(completed with lock and handle).
- Dimension** : 1200 (L) x 450 (W) x 1700 (H) mm.
- Castors** : Four 3" castors, two with brakes.
- Finishing** : Powder coating.

Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com
www.grahateknomedika.com

MADE IN INDONESIA

