



INSTRUMENT CABINET

1 DOOR
8211401



SPECIFICATIONS :

- | | |
|---------------------|---|
| Construction | : Main frame made of steel pipe.
Upper frame made of steel plate.
Three shelves made of glass (5 mm).
Upper door made of square steel pipe,
(completed with lock and handle).
Lower door made of steel plate,
(completed with lock and handle). |
| Dimension | : 700 (L) x 450 (W) x 1700 (H) mm. |
| Castors | : Four 3" castors, two with brakes. |
| Finishing | : Powder coating. |

Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com

www.grahateknomedika.com

MADE IN INDONESIA

