



STRETCHER BACK RAISE

8211002

SPECIFICATIONS :

- Construction** : Upper frame made of steel pipe.
 Back rest frame made of steel pipe.
- Dimension** : 2000 (L) x 600 (W) x 800 (H) mm.
- Adjustment Range** : Back rest adjustable by gas spring up to 70°.
- Castors** : Four 4" castors, two with brakes.
- Finishing** : Polished.
- Accessories** : IV stand made of stainless steel.
 Mattress made of foam covered by vinyl.
 Multi-purpose tray.



Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
 Jakarta Selatan 12240 - Indonesia
 P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
 Sawangan, Depok - Indonesia
 P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com

www.grahateknomedika.com

MADE IN INDONESIA

