



VERLOS BED

STAINLESS STEEL

8221101



SPECIFICATIONS :

Construction

- : Main frame made of stainless steel pipe.
- Back rest frame made of stainless steel pipe.
- Extension frame made of stainless steel pipe.

Dimension

- : Overall 2100 (L) x 700 (W) x 1000 (H) mm.

Adjustment Range

- : Adjustable by gas spring up to 70°.

Castor

- : Four 4" castors, two with brakes.

Finishing

- : Polished.

Accessories

- : Mattress made of foam covered by vinyl.
- Gynaecology leg strut.
- IV stand made of stainless steel.
- Bowl made of stainless steel.



Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com

www.grahateknomedika.com

MADE IN INDONESIA

