



X - RAY VIEWER

1 SCREEN

8211901



SPECIFICATIONS :

Construction :

Main frame made of steel plate and acrylic.

Dimension :

400 (L) x 125 (W) x 550 (H) mm.

Power Supply :

220 V / 50 Hz.

Lamp :

1 x 32 Watt (On/Off switch).

Finishing :

Powder coating.

Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com

www.grahateknomedika.com

MADE IN INDONESIA

