



GYNAECOLOG TABLE 8211009

SPECIFICATIONS :

- | | |
|-------------------------|---|
| Construction | : Mattress base consists of 3 parts of the steel box frame. |
| Dimension | : 200 (L) x 60 (W) x 80 (H) cm. |
| Adjustment Range | : Each position can be set manually. |
| Castor | : Four 4" castors, two with brakes. |
| Finishing | : Powder coating. |
| Accessories | : Mattress made of foam covered by vinyl.
Gynaecology leg strut.
IV stand made of stainless steel.
Bowl made of stainless steel. |



Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com

www.grahateknomedika.com

MADE IN INDONESIA

