



**GRAHA TEKNO MEDIKA**

Hospital Furniture & Equipment Manufacturing



## LED X-RAY VIEWER 1 SCREEN AB1911



### SPECIFICATIONS :

- Construction** : Main frame made of steel plate, aluminium and acrylic.
- Lamp** : LED 24 VDC.  
Power consumption 20 watt.  
On/Off push button.  
Adjustable light dimmer.
- Power supply** : 220V / 50 HZ.
- Brighness** : 10 - 3500 lux.
- Dimension** : 535 (L) x 495 (H) x 30 (W) mm
- Visible Size** : 370 (L) x 422 (H) mm
- Finishing** : Powder coating.



#### Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24  
Jakarta Selatan 12240 - Indonesia  
P. +6221 7264150 | F. +6221 7264152

#### Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,  
Sawangan, Depok - Indonesia  
P. +6221 7494512 | F. +6221 7494512

[info@grahateknomedika.com](mailto:info@grahateknomedika.com)

[www.grahateknomedika.com](http://www.grahateknomedika.com)

**MADE IN INDONESIA**

