



GRAHA TEKNO MEDIKA

Hospital Furniture & Equipment Manufacturing



INFUSION PUMP APF3001



SPECIFICATIONS

Overall Dimension (mm)	: 260 x 130 x 91 mm
Weight	: 1.8 kg
Pumping Mechanism	: Stepping motor drives the drive shaft.
Flow rate	: 0.01 - 1500 ml/h
Accuracy	: $\pm 3\%$ (in room temperature 25° C, with water)
Range Total Volume	: 0.01 - 9999 ml
Range Infused Volume	: 0.01 - 9999 ml
Bolus range	: 1 - 1500 ml/h (default 700 ml/h)
Bolus volume	: 0.01 - 99.9 ml
Occlusion Pressure	: 60 - 850 mmHg
Display	: Mono Graphic LCD, (resolution 240 x 64)

Power :

AC 100 - 240 VAC, 50 - 60 Hz

DC 12 Volt

Ni - MH Rechargeable

Battery 2000 mAh

Consumption 34 VA

Protection :

Class I, Electric Shock.

Type CF, Capacity for Electric.

IPX2, Liquid Infiltration.

Alarm & Sensor :

Air Bubble Alarm, Occlusion Alarm & Sensor,
Battery Alarm, Power Alarm, Near Complete Alarm,
Nurse Call, Air Sensor, Drop Sensor.

Marketing Office :

Jl. Sultan Iskandar Muda Kav. 24
Jakarta Selatan 12240 - Indonesia
P. +6221 7264150 | F. +6221 7264152

Factory & Main Logistic :

Jl. Abdul Wahab No. 59, Cinangka,
Sawangan, Depok - Indonesia
P. +6221 7494512 | F. +6221 7494512

info@grahateknomedika.com

www.grahateknomedika.com

MADE IN INDONESIA

